



Fig.1 H301-HM-ACTIVE



Fig.2 H301-HM-ACTIVE + CO2 UNIT BL

Technical Specifications

Relative Humidity range	50% - 95 %
Humidity sensor resolution	1 %
Humidity sensor accuracy	±2 %
Humidification technique	Water permeable membrane immersed in water at controlled temperature
Dry gas input	From any Bold Line Gas Controller
Humid gas output	In temperature controlled tube
Compatible with	CO ₂ O ₂ UNIT BL [0-10; 1-18] CO ₂ O ₂ UNIT BL [0-20; 1-95] CO ₂ O ₂ UNIT BL [0-10; 0-1]
Power	55 W
Voltage	24 V
Certification	CE

Features

- Feedback From Humidity Sensor
- Control of Relative Humidity inside Chamber
- Water Permeable Membrane
- Heated connection tube to prevent condensation
- Water Heater
- OKO Touch operated




Fig.3 Humidity Controller OKO Touch Page

Assessment of Humidification Performance

Methods

Humidification performance was assessed by measuring water evaporation in a Top Sage Incubator loaded with 4 35mm Petri Dish (with no lid), each one initially containing 3 ml of water. The top stage incubator was kept at 37°C and was positioned on a 1 mg precision balance, connected to a PC for data acquisition. A continuous flux of 0.4 l/min of 5% CO₂ in air was sent to the Top Stage Chamber after being humidified by the Active Humidity Controller. Relative Humidity inside the chamber was varied from 55% to 95%.

Experimental apparatus

Figure 4 reports a picture of the experimental apparatus. The equipment tested was (numbers in parenthesis refer to Figure 4):

- Top Stage Chamber: H301-EC-BL (1); T set at 37.0°C
- Specimen Holder: H301-EC-4x35
- Petri Dish: 4 x 35 mm Petri Dish
- Initial water content: 3 ml in each Petri Dish
- T controller: H301-T-UNIT-BL (3)
- CO₂ Controller: CO2 UNIT BL (3)
- Gas Flow rate: 0.4 l/min, at 5.0% CO₂ in air
- Humidity Controller: H301-HM-ACTIVE (4)
- Balance: Precisa Gravimetrics AG-XB 620 M (2)



Fig.4 Experimental apparatus

Results

Figure 5 reports the weight loss of the Top Stage Chamber due to evaporation as a function of time. Relative Humidity was varied from 55% to 95%. Data were shifted to the same origin for the sake of clarity. The Dry Gas data were obtained by by-passing the Active Humidity Controller.

Evaporation rate (\hat{e}) reported in Figure 5 was calculated by dividing the weight loss (in micrograms per hour) by the air /water contact area (surface of four petri dish). Data are therefore reported as $\mu\text{g}/\text{hr}/\text{mm}^2$.

Data show that evaporation in the chamber becomes negligible when Relative Humidity in the chamber is set to 95%.




Fig.5 Weigh loss due to evaporation in a Top Stage Incubator loaded with 4 35mm Petri Dish as a function of Relative Humidity inside the chamber